

# [1002] Correlation of nano-scale electrical and topographical mapping of buried nanoscale semiconductor junctions

Ghazi Alsharif [UK]<sup>1</sup>, Linda Hanel [Germany]<sup>2</sup>, Alex Robson [UK]<sup>1</sup>, Jorg Schulze [Germany]<sup>2</sup>, Benjamin Robinson [UK]<sup>1</sup>, Oleg Kolosov [UK]<sup>1</sup>

Lancaster University<sup>1</sup>, Institute of Semiconductor Engineering, University of Stuttgart<sup>2</sup>


In recent years, germanium tin (Ge<sub>0.9</sub>Sn<sub>0.1</sub>) group-IV semiconductor materials attracted interest as promising candidates for inclusion in optoelectronic devices. They provide a possible route for realising direct-bandgap semiconductors with increasing Sn content [1], potentially allowing to create LED in a non iii-v semiconductor materials. With the continuing demand for miniaturisation in modern electronic devices, there is a linked demand for new characterisation techniques operating in the nanoscale regime.

Here we show that by combining an Ar-ion cross-sectioning technique and a customised scanning probe microscopy (SPM), one can extract and explore the electrical properties of the subsurface Ge<sub>0.9</sub>Sn<sub>0.1</sub> structures with a resolution of ~30 nm. This high spatial resolution current mapping is needed to fully understand the nanoscale mechanisms of electrical transport these complexly structured semiconductor nanostructures and assist in the development of new electronic nano-devices.

Beam-Exit X-section Polishing (BEXP) [2], Fig. 1 (a), was used to form high quality cross sections of the layered Ge<sub>0.9</sub>Sn<sub>0.1</sub> samples enabling direct access with scanning spreading resistance microscopy (SSRM) electrical mapping. Our SSRM is based on a conventional SPM system (Bruker Multimode with Nanoscope 4a), with a highly-electrically conductive diamond probe and custom probe-signal pre-amplifier and signal access capabilities. The setup enabled the measurements and interpretation of the nanoscale electrical transport properties and mapping of spreading resistance of the subsurface layers under ambient environment conditions.

In this paper we report nanoscale measurements on 100 and 200 nm layer of Ge<sub>0.9</sub>Sn<sub>0.1</sub> on the gradient Ge virtual substrates (GeVS) MBE grown on the Si wafer. We obtained SSRM images, as well as captured (I-V) curves with and without illumination.

The annealing of the GeVS allowed low density of the defects immediately under the top GeSn<sub>x</sub> layer. SSRM allowed to observe that subsequent annealing step of the total structure provide notable variation of the conductivity of the GeSn<sub>x</sub> layer, as well as to uncover the high defect density layer of the GeVS immediately above Si substrate. These were confirmed by the concurrent nanomechanical mapping via ultrasonic force microscopy (UFM). In conclusion, combination of precise sectioning via BEXP and subsequent SSRM imaging that was demonstrated for the first time in this paper, will be a powerful tool for the investigation of the multilayer semiconductor heterostructures.



([https://www5.shocklogic.com/Client\\_Data/RMS/al/MMC2017/upload/IRMS-MMC2017-37077841-Fig.1.jpg](https://www5.shocklogic.com/Client_Data/RMS/al/MMC2017/upload/IRMS-MMC2017-37077841-Fig.1.jpg))

Fig 1. (a) BEXP cross section which produces a shallow angled ( $5 - 15^\circ$ ) cut through the sample, and enables subsurface SPM measurements [3], (b) scanning spreading resistance microscopy (SSRM) image of the tip-material conductance, (c) composite image with fixed AC voltage (1Vpp) with different DC offset as indicated on the left of the image, (d) I-V curves taken at points in (b), (e)  $dI/dV$  of (d).